


Centrolock® – System Tray Layout

centroLock®
Guided Transverse Osteotomy System



Centrolock Instruments

	Part No	Description	Qty
A	XMS01033	Retractor	2
B	ACC1001P0020	Centrolock K-wires holder	1
C	CKW01012	K-wire Ø2 lg. 100 TR-RD extra sharp*	4
D	CKW01013	K-wire Ø2 lg. 150 TR-RD extra sharp*	4
E	XPP01005D	Centrolock temporary fixation pin*	2
F	XMS01028	Centrolock cutting guide	1
G	XMS01029	Centrolock elevator	1
H	XMS01009	Percutaneous rasps	Optional
I	XRE01014	Centrolock cannulated reamer	1
J	XTI06012	Centrolock trial implant – Left step 2mm	1
K	XTI06014	Centrolock trial implant – Left step 4mm	1
L	XTI06016	Centrolock trial implant – Left step 6mm	1
M	XTI06022	Centrolock trial implant – Right step 2mm	1
N	XTI06024	Centrolock trial implant – Right step 4mm	1
O	XTI06026	Centrolock trial implant – Right step 6mm	1
P	XMS01030	Centrolock impactor	1
Q	XHA01001	AO handle	1
R	XHA01002	AO Ratcheting Handle	1
S	XDG01019	Centrolock locking drill guide – screw Ø2.5mm	2
T	XDB01020D	Centrolock drill bit Ø1.8mm*	2
U	XGA01011	Centrolock depth gauge	1
V	XSD01003	Centrolock screwdriver tip T7	2
W	XMS01026	Centrolock targeting / compressing guide	1
X	XGA01003	Screw measurer	1
Y	XDG01018	Centrolock drill guide – screw Ø2mm	2
Z	XDB01019D	Centrolock drill bit Ø1.5mm*	2
AA	XRE01022	Centrolock Countersink	1
-	XMS01036	Centrolock straight impactor	Optional